



# NoDrip™ Sterile Wrap & Instrument Protection




## NoDrip™

### Sterile Wrap Protection, Steam Sterilization Compatible


Use NoDrip Padding between the tray and the sterile wrapping to prevent rips and tears caused by a tray's sharp edges or protruding feet. The absorbent lint-free paper lining helps wick moisture away from the CSR wrap while the foam padding provides cushioning.

### Protect delicate instruments with NoDrip Padding

It can be used to cushion and protect delicate instruments in a tray. Ideal for heavy instruments sets.



SKU	DESCRIPTION	SIZE	UNITS/CASE
ND4500	NoDrip Roll - Easy Tear 3" Perforation	15" x 40'	2 Rolls
ND4550	NoDrip Roll - Easy Tear 3" Perforation	15" x 100'	2 Rolls
ND4600	NoDrip Roll - Easy Tear 3" Perforation	12" x 32.5'	4 Rolls
ND4650	NoDrip Roll - Easy Tear 3" Perforation	9" x 32.5'	4 Rolls
ND124	NoDrip End Caps	12" x 4"	150 End Caps



### NoDrip Acrylic Roll Dispenser

Order#: NDRD15M  
Size: 16"W x 19.5"D x 11.25"H

### Protect Sterile Wrapping with Steam Sterilization Compatible Tray Corner Guards

Tray Corner Guards are constructed of urethane foam and act as padding for the tray. Compatible with steam sterilization.

Order#: SC8403  
Size: 3" x 3"  
Units/Case: 288 Guards

For sterile wrap protection and tray padding compatible with low-temperature plasma sterilization see Plasma-Cel™ and BeeSafe™ products on pages 7-8.



Cygnus Medical 800.990.7489 | [www.cygnusmedical.com](http://www.cygnusmedical.com)

965 West Main Street, Branford, CT 06405 | Email: [sales@cygnusmedical.com](mailto:sales@cygnusmedical.com)