



Soiled instruments leaving a procedure must be in a closed container labeled as biohazard.



## SingleCycle™ Instrument Transport Trays and Biohazard Labels

SingleCycle provides secure, tamper-proof protection for soiled instruments. The single-use containers are designed to protect even the most delicate instruments from damage, and are available in a variety of sizes. Not only do they eliminate the need for manual disinfection of reusable bins, they can help to reduce cross-contamination in the work area.

**Meets the SGNA, AAMI & AORN standard of a closed, puncture-resistant container that protects staff and patients from potentially infectious organisms.**

- Made of a 100% renewable biodegradable plant-based material.
- Single-use, space efficient, stackable, and nestable.
- Autoclavable, can be steam sterilized to transport sterile wrapped trays.

[See side two for tray options and ordering information >](#)



**Cygnus Medical 800.990.7489 | [www.cygnusmedical.com](http://www.cygnusmedical.com)**

965 West Main Street, Branford, CT 06405 | Email: [sales@cygnusmedical.com](mailto:sales@cygnusmedical.com)



**--- SINGLE-USE BIODEGRADABLE AUTOCLAVABLE ---**

**SingleCycle™ Disposable Instrument Trays and Biohazard Labels make it easy to properly transport & label wrapped instruments traveling to clinics.**

SKU	DESCRIPTION	SIZE	UNITS/CASE
IT0308	Small SingleCycle™ Instrument Transport Container	Inner Dimension 3" x 8" Outer Dimension 4.4" x 9.45"	100 Trays
IT0308-SI	Bundle - small SingleCycle™ Instrument Transport Container (IT0308) & Biohazard Label (adhesive backing)	Same dimensions as IT0308 and SCB15	100 Trays 100 Labels
IT0612	SingleCycle™ Instrument Transport Container	Inner Dimension 6" x 12" Outer Dimension 8" x 14"	100 Trays
IT0612-SI	Bundle - SingleCycle™ Instrument Transport Container (IT0612) & Biohazard Label (adhesive backing)	Same dimensions as IT0612 and SCB20	100 Trays 100 Labels
IT0713	SingleCycle™ Instrument Transport Container	Inner Tray Dimensions 7" x 13" Outer Tray Dimensions 10" x 15.5"	60 Trays
IT0713-SI	Bundle - SingleCycle™ Instrument Transport Container (IT0713) & Biohazard Label (adhesive backing)	Same dimensions as IT0713 and SCB25	60 Trays 60 Labels
SC0608	Foam Tray Liners for small SingleCycle™ (IT0308) Instrument Transport Container.	6" x 8"	100 Liners
SC1212	Foam Tray Liners for SingleCycle™ (IT0612) Instrument Transport Container.	12" x 12"	100 Liners
SCB15	Biohazard label, used with IT0308 SingleCycle Trays	15" x 1.75"	500 labels
SCB20	Biohazard label, used with IT0612 SingleCycle Trays	20" x 1.75"	500 labels
SCB25	Biohazard label, used with IT0713 SingleCycle Trays	25" x 1.75"	500 labels



**Cygnus Medical 800.990.7489 | www.cygnusmedical.com**

965 West Main Street, Branford, CT 06405 | Email: sales@cygnusmedical.com