


# Disposable Endoscopic Channel Brushes

## Single-Use and disposable channel cleaning brushes.

Cygnus Medical offers a variety of endoscopic channel cleaning brushes designed to meet your endoscope reprocessing demands. Nylon bristles clean efficiently and are safe for use in endoscope channels.

### CB1015



### CB1015

#### Disposable Endoscopic Channel Brush

Compatible Channel Sizes: 1.0 mm - 1.5 mm


Brush Diameter: 2 mm

Brush Length: 10 mm

Total Length: 1000 mm

Units/Case: 100 Brushes

### CB2024



### CB2024

#### Disposable Endoscopic Channel Brush

Compatible Channel Sizes: 2.0 mm - 2.4 mm


Brush Diameter: 3 mm

Brush Length: 20 mm

Total Length: 1800 mm

Units/Case: 50 Brushes

### CB2024-L



### CB2024-L

#### Disposable Endoscopic Channel Brush

Compatible Channel Sizes: 2.0 mm - 2.4 mm


Brush Diameter: 3 mm

Brush Length: 20 mm

Total Length: 2300 mm

Units/Case: 50 Brushes

### CB2845



### CB2845

#### Disposable Endoscopic Channel Brush

Compatible Channel Sizes: 2.8 mm - 4.5 mm


Brush Diameter: 5 mm

Brush Length: 20 mm

Total Length: 2300 mm

Units/Case: 50 Brushes

### VBD1005



### VBD1005

#### Disposable Double-Ended Valve Brush


Brush Diameters: 10 mm and 5 mm

Brush Lengths: 30 mm and 20 mm

Total Length: 150 mm

Units/Case: 100 Brushes

### BB10



### BB10

#### Bristle Button Brush

Brush Diameter: 3 mm

Brush Length: 25 mm

Total Length: 102 mm

Units/Case: 100 Brushes

(25 Units per Bag, 4 Bags per Case)

All brush diagrams are for information only. Not drawn to scale.

