



# PathoGene™ qPCR Detection Kits

Primers  
Probes  
Master Mix

**All you need in one little bead.**  
qPCR prep time in under 5 minutes.

## PathoGene™ qPCR Detection Kits

<b>C. auris Detection</b> PDK-CAU08-T	<b>C. diff Detection</b> PDK-CDF08-T	<b>Norovirus Detection</b> PDK-NOR08-T	<b>MRSA Detection</b> PDK-MRS08-T
--	---	---	--------------------------------------



**Complete pre-plated PCR mix in a single freeze-dried bead that includes primers, probes and master mix.**

- PathoGene Detection Kits streamline sample preparation, while removing the potential for human error.
- High analytical sensitivity and specificity. Returns consistent, accurate results.
- PathoGene™ Transport Media inactivates multidrug-resistant pathogens at point-of-care, eliminating the need for an RNA extraction step.
- Eliminates cold chain storage. Specimen sample collection, storage, and transport can be conducted at room temperature, with the samples remaining stable for several weeks.

See the complete line of PathoGene Products & ordering information on the next page.

**RUO** Research Use Only. Not for use in Diagnostic Procedures.

SKU	PRODUCT	CONTENTS
<b>PathoGene qPCR Detection Kits</b>		
PDK-CAU08-T	C. auris* qPCR Detection Kits	8 Oligo Reaction Beads (1 Strip) 8 Vials Pathogene Transport Media 8 Swabs
PDK-CDF08-T	C. diff** qPCR Detection Kits	8 Oligo Reaction Beads (1 Strip) 8 Vials Pathogene Transport Media 8 Swabs
PDK-NOR08-T	Norovirus qPCR Detection Kits	8 Oligo Reaction Beads (1 Strip) 8 Vials Pathogene Transport Media 8 Swabs
PDK-MRS08-T	MRSA*** qPCR Detection Kits	8 Oligo Reaction Beads (1 Strip) 8 Vials Pathogene Transport Media 8 Swabs
<b>PathoGene Test Strips</b>		
PDK-CAU08	C. auris* qPCR Detection Test Strip	1 Strip/8 Reactions
PDK-CDF08	C. diff** qPCR Detection Test Strip	1 Strip/8 Reactions
PDK-NOR08	Norovirus qPCR Detection Test Strip	1 Strip/8 Reactions
PDK-MRS08	MRSA*** qPCR Detection Test Strip	1 Strip/8 Reactions
<b>Positive/Negative Controls</b>		
PPNC01-12	Positive/Negative Controls - for use with Pathogene C. auris, C. diff, Norovirus, and Detection Strips.	1 Strip/8 Reactions
<b>PathoGene Transport Media</b>		
PTM-50	Pathogene Transport Media	1 Box of 50 Vials

\*Candida auris (C. auris)

\*\*Clostridioides difficile (C. diff)

\*\*\*Methicillin-resistant Staphylococcus aureus (MRSA)