


# NoDrip™ Padding



NoDrip Full Sheet Protection



NoDrip End Caps



NoDrip In Tray



NoDrip Roll Dispenser

## NoDrip™ Protection

NoDrip Padding can be used to protect delicate instruments in a tray or be used beneath a tray to prevent rips and tears caused by a tray's sharp edges or protruding feet. The absorbent lint-free paper lining helps wick moisture away from the instruments or CSR wrap while the foam padding provides cushioning. Ideal for heavy instruments sets.

## NoDrip™ End Caps

End Caps focus on the "high heel" effect of pressure points caused by the tray's feet or sharp edges. One or more layers can be used to prevent rips and tears in the CSR wrap.

*Call Cygnus Medical today or visit our website to learn more about NoDrip Products.*

SKU	DESCRIPTION	SIZE	UNITS/CASE
ND4500	NoDrip Roll – Easy Tear 3" Perforation	15" x 40'	2 Rolls
ND4550	NoDrip Roll – Easy Tear 3" Perforation	15" x 100'	2 Rolls
ND4600	NoDrip Roll – Easy Tear 3" Perforation	12" x 32.5'	4 Rolls
ND4650	NoDrip Roll – Easy Tear 3" Perforation	9" x 32.5'	4 Rolls
ND124	NoDrip End Caps	12" x 4	150
NDRD15	NoDrip Roll Dispenser	16"w x 18"d x 7"h	1



**CYGNUS**  
M E D I C A L

1.800.990.7489  
sales@cygnusmedical.com  
www.cygnusmedical.com