



Scope Transport and Protection

Single-Use Rigid Containment



Lapses in the cleaning of reusable endoscope transport trays can lead to cross-contamination in the surgical work area. Bedside cleaning is performed within the Oasis Tray, keeping residual run-off clearly differentiated clean and soiled scopes. The Oasis Scope Transport Tray offers an economical single-use alternative that protects both patients and staff from contamination and infection.

Oasis and Oasis PLUS Trays are Latex-Free.



OT1230



OT1016



OP1330

Oasis® Scope Transport Trays

SKU	DESCRIPTION	SIZE	UNITS/CASE
OT1230	Oasis Biodegradable Tray with Reversible OR Lid.	20.5" x 16.5"	30 Trays 30 Lids
OT1016	Oasis Biodegradable Tray with Reversible OR Lid. Perfect for small diameter scopes and small instruments.	16" x 10"	50 Trays 50 Lids
OP1330	Oasis PLUS Biodegradable Tray with Reversible OR Lid.	21.5" x 19.5"	30 Trays 30 Lids

Single-Use and disposable.

Made of a biodegradable plant-based material.

Space efficient, stackable and nestable.

Dispensing, Collection and Transport – Side 2

Oasis® and Oasis® PLUS Dispensing and Collection



Oasis and Oasis PLUS Dispensing Dolly

Load clean nested trays in the Oasis Dispensing Dolly. The cart features a stainless steel base, 50 tray capacity with lid shelf, removable prep table top, non-marking swivel casters. Dolly color - green.



Oasis and Oasis PLUS Collection Unit

Safely dispose of soiled trays in the Oasis Collection Unit. The unit features a stainless steel base, non-marking swivel casters. Unit color - red.

Biodegradable Liners

Earth-Friendly liners sized for the Oasis Collection Unit.



SKU	DESCRIPTION	SIZE	UNITS/CASE
DD40-G	Oasis/Oasis PLUS Dispensing Dolly, color-green.	23"W x 21"D x 36"H	1

SKU	DESCRIPTION	SIZE	UNITS/CASE
CU10-R	Oasis/Oasis PLUS Collection Unit, color - red.	23"W x 21"D x 36"H	1
4253-100	Biodegradable Liners	Fits CU10-R	100 Liners

Oasis® and Oasis® PLUS Transport Carts

All transport carts are available in both **green/Clean** and **red/Soiled**.



TC04-G



TC06-R

Oasis and Oasis PLUS Transport Carts are available 4, 6, or 8 tray capacities.

Oasis Trays fit perfectly in the Oasis Transport Cart for secure and easy transport. All carts are available in green (clean scope cart) and red (soiled scope cart) with non-marking swivel caster wheels.

SKU	DESCRIPTION	SIZE	UNITS/CASE
TC04-G	4 tray Oasis/Oasis PLUS Transport Cart, color - green.	23"W x 21"D x 32"H	1
TC04-R	4 tray Oasis/Oasis PLUS Transport Cart, color - red.	23"W x 21"D x 32"H	1
TC06-G	6 tray Oasis/Oasis PLUS Transport Cart, color - green.	23"W x 21"D x 44"H	1
TC06-R	6 tray Oasis/Oasis PLUS Transport Cart, color - red.	23"W x 21"D x 44"H	1
TC08-G	8 tray Oasis/Oasis PLUS Transport Cart, color - green.	23"W x 21"D x 56"H	1
TC08-R	8 tray Oasis/Oasis PLUS Transport Cart, color - red.	23"W x 21"D x 56"H	1
TP-50	Transport Cart Top Cover	23"W x 21"D	1

