



Ask your Sales Rep about Bundle Pricing. Free Brush Organizer Included!

Pushing



Pulling

The pull-through design increases the amount of brushing action without causing buckling.



pulling vs. pushing

Cygnus Medical Pull-Through Channel Brushes

Our nylon Pull-Through Channel Brush is designed to significantly increase the amount of brushing action without causing buckling. The flexible leader makes it easy to clean hard-to-reach narrow channels. Simply insert the leader into the channel and pull the brush through. Use with Simple2™ Enzymatic Detergent for the most effective pre-clean process.

The Pull-Through Channel Brush can be used on all laparoscopic and cannulated stainless steel instruments.

SKU	BRUSH SIZES	LENGTH	UNITS/PACK
15-010	1 mm (3 French)	15" Leader	5 Brushes
15-013	1.33 mm (4 French)	15" Leader	5 Brushes
15-016	1.6 mm (5 French)	15" Leader	5 Brushes
15-020	2 mm (6 French)	15" Leader	5 Brushes
16-030	3 mm (9 French)	16" Leader	5 Brushes
16-033	3.3 mm (10 French)	16" Leader	5 Brushes
16-040	4 mm (12 French)	16" Leader	5 Brushes
M16-400	4 mm (12 Fr.) Coarse Bristles	16" Leader	5 Brushes
M16-600	6 mm (18 Fr.) Soft Bristles	16" Leader	5 Brushes
CB406-PT	Pull-Through Brush Organizer	Call for Size	1 Organizer