


Tray beltTM Sterile Wrap Protection



Compatible with both steam and low-temperature sterilization.




Available In rolls or as pre-cut belts.



The Tray Belts Dispenser provides easy access to the belts when you need to apply them to trays.



Belt application outside of sterile wrapping cannot only help to prevent external damage that can occur during sterilization, storage, and transport, but may also prevent moisture problems like wet packs.



SKU	DESCRIPTION	SIZE	UNITS/CASE
Tray Belts Compatible with Steam Sterilization			
TBR3100NP	Tray Belts Roll	3" X 100'	8 Rolls and 1 Roll of Indicator Tape.
TBPC30	30" Pre-Cut Belts in a Dispensing Bag.	3" X 30"	200 Belts (4 Bags of 50 Belts per Case.)
TBPC46	46" Pre-Cut Belts in a Dispensing Bag.	3" X 46"	200 Belts (4 Bags of 50 Belts per Case.)
Tray Belts Compatible with Low-Temperature Sterilization			
HCB0392	Honeycomb Tray Belts	3" X 92'	4 Rolls and 1 Roll of Indicator Tape.
TBLT-RL	Tray Belts for Low-Temperature Sterilization (Sterrads and V-Pro)	3" X 100'	2 Rolls and 1 Roll of Indicator Tape.
Dispensers			
TBRD10	Tray Belts Acrylic Roll Dispenser	7.5"w X 20"d X 11"h	1 Dispenser
TBRDMETAL	Tray Belts Metal Roll Dispenser	7.5"w X 20"d X 11"h	1 Dispenser



Cygnus Medical
 965 West Main Street, Branford, CT 06405
www.cygnusmedical.com | 800.990.7489
 Patent Number 10751436



**Learn More
 Watch the
 Video.**