



## The information you need at a glance.



Easily attaches to itself with pressure-sensitive adhesive.



SGNA supports a 7-day storage interval for reprocessed endoscopes, but only if they were reprocessed and stored according to professional guidelines and manufacturer instructions.

**Clearly labeled scope tags help your staff stay in compliance. Quickly identify scopes that require reprocessing with easy-to-read color-coded tags.**

### DAILY COUNTDOWN™ Scope Tags

Pressure-sensitive adhesive strip on the tag back.  
12" x 1.75" | 400 Tags per Case (4 Packs of 100)

WEEKDAY	ORDER #
Monday	CDM610
Tuesday	CDT620
Wednesday	CDW630
Thursday	CDTH640
Friday	CDF650
Saturday	CDS660
Sunday	CDSUN670

#### Combo Pack ORDER # CDCOMB700

700 Tags per Case  
(7 Packs of 100,  
1 Pack of Each Day)

### COMPLIANCE STANDARDS

## ANSI/AAMI ST91:2015

### Flexible and Semi-Rigid Endoscope Reprocessing in Health Care Facilities

#### 10. Storage of Reprocessed Endoscopes

##### 10.1 General Considerations

“Attach a tag or label (or other method) to document that the scope has been cleaned or high-level disinfected. The tag or label is affixed to the endoscope after it has been processed and before it is placed in the storage cabinet.”

Contact Cygnus Medical at 800.990.7489 to learn more about DAILY COUNTDOWN Scope Tags.