

Hoping instrument jaws stay open during sterilization cycles may not be enough.

Achieve 🤞 of mind with the Gator Instrument Clip.

Exposing instrument surface area during sterilization is critical in order to ensure sterility. Gator Instrument Clips from Cygnus Medical help to keep instrument jaws open during the sterilization cycle.

Peel pouch instruments must be sterilized on their side edges. In this position gravity pulls the instrument jaws closed. Gator Instrument Clips prevent this from happening.

SKU	DESCRIPTION	SIZE	PER CASE
GIC190-5	Use with Ring-Handled Scissors	1.9"	1 Bag of 500
GIC190-10	Use with Ring-Handled Scissors	1.9"	1 Bag of 1000
GIC270-5	Use with Ratcheted Ring-Handled Instruments	2.7"	1 Bag of 500
GIC270-10	Use with Ratcheted Ring-Handled Instruments	2.7"	1 Bag of 1000




- Latex-Free
- Single-Use
- Multiple Sizes



Cygnus Medical

965 West Main Street, Branford, CT 06405
800.990.7489 www.cygnusmedical.com



The Gator Instrument Clip™
PATENT PENDING