



Wish you could reduce endoscope channel drying times?

Now you can with Airtime™ Drying Swabs



Airtime Channel Swabs

Airtime Channel Swabs virtually remove all moisture inside the biopsy/suction channel in one pass. Constructed with multiple foam heads, this novel approach drastically reduces drying times for flexible endoscopes.

The leading non-absorbent head effectively expels excess water, while the trailing absorbent foam head captures any remaining moisture.



Airtime Port Swabs

Airtime Port Swabs remove water trapped and pooled inside valve chambers.

Drying Swab Ordering Information

SKU	Description	Size	Units/Case
ADS03-100	Channel Swab, non-sterile	3mm* OD	100
ADS04-100	Channel Swab, non-sterile	4mm* OD	100
ADS05-100	Channel Swab, non-sterile	5mm* OD	100
APS-200	Port Swab, non-sterile	N/A	200

*nominal



CONTACT CYGNUS MEDICAL

965 West Main Street, Branford, CT 06405
Phone: 800.990.7489 Email: sales@cygnusmedical.com
www.cygnusmedical.com

This product is patent-pending