



clear path™

Polyester Wipers for Drying

High-Absorbency Wipers

Drying flexible endoscopes after high-level disinfection is a required and critical step in the scope reprocessing cycle. Clear Path high-absorbency lint-free wipers are available in sterile and non-sterile packaging, and provide fast, effective removal of moisture without transferring lint onto endoscopes.

FEATURES

- Highly absorbent polyester.
- Sterile and non-sterile options.
- Lint and particulate-free.
- Various packaging configurations available.

APPLICATIONS

Use with the Airtime Instrument Channel Dryer. Dry the endoscope exterior with absorbent wipers and use Airtime to dry scope channels.

SKU	DESCRIPTION	SIZE	UNITS/CASE
ST0909	Clear Path™ lint-free, absorbent, STERILE wipers, 2 wipers per bag - 100 bags per case, 200 wipers total. Polyester, latex-free, and single-use.	9" x 9"	200
WN0909-2	Clear Path™ lint-free, absorbent, NON-STERILE wipers, 2 wipers per bag - 100 bags per case, 200 wipers total. Polyester, latex-free, and single-use.	9" x 9"	200
PW0909100	Clear Path™ lint-free, absorbent, NON-STERILE wipers, 100 wipers per bag - 6 bags per case, 600 wipers total. Polyester, latex-free, and single-use.	9" x 9"	600

Clear Path Microfiber Wipers for Cleaning - Side 2 >





Microfiber Wipers for Cleaning

STERILIZATION CONTAINER DE-LINTING

The dust and lint that settles on lids and bases of closed sterilization containers can migrate onto surgical instruments during steam sterilization. Clear Path Wipers remove dust and particulates prior to sterilization, providing a higher level of patient care.

SKU	DESCRIPTION	SIZE	UNITS/CASE
FS4193	Disposable wipers with heat-sealed edges, constructed of woven, lint-free, particulate-free microfiber, 25 wipers per bag - 4 bags per case.	11" x 11"	100

SURFACE CLEANING

Clear Path split microfiber technology captures spores and other viable pathogens at a microscopic level. These single-use wipers prevent the cross-contamination of surface areas, and can be used with any brand of surface disinfectant. Ideal for use in operating rooms, intensive and cardiac care units, neonatal units, sterile processing departments and procedure rooms.

SKU	DESCRIPTION	SIZE	UNITS/CASE
W5412 - 30	Disposable wipers, constructed of non-woven, lint-free, low-particulate microfiber, 30 wipers per waterproof dispenser package - 20 packages per case.	8" x 11.75"	600
W5412 - 150	Disposable wipers, constructed of non-woven, lint-free, low-particulate microfiber, 150 wipers per bag - 8 bags per case.	10" x 12"	1200

