

Protect your probe, from itself.

Prevent cross-contamination with Dual-Cavity Transport Bags.



OPEN AREA (not sealed)
Allows for passage of the
instrument cords.

HEAT-SEALED STRIP
Separates the **BIOHAZARD**
compartment from the
CORD compartment.

Order#: DCB2017
Size: 17" x 20"
Case Quantity: 100 Bags

Not all parts of internal probes require high-level disinfection (HLD).

Cords and connectors need to be kept separate from the internal transducer that has been exposed to gross contamination.

DuPlex™ Dual-Cavity Instrument Transport Bags isolate dirty portions of instruments and prevent the cross-contamination of sections that are not high-level disinfected.

DuPlex™

DUAL-CAVITY INSTRUMENT TRANSPORT BAG

